

capillus^{rx}

BY PHYSICIANS ONLY

HELPS PREVENT HAIR LOSS
HELPS HAIR REGROWTH



MOBILE LASER THERAPY CAP



Laser Density

1560

mW Total Output

5

Year Warranty

Maximum Coverage Physician-Only

- ✓ REVERSES THE PROCESS OF HAIR LOSS
- ✓ SAFE & EFFECTIVE FOR TREATING HAIR LOSS DUE TO HEREDITY
- ✓ TREATS THINNING HAIR AND BALDING
- ✓ NURTURES FULLER, HEALTHIER LOOKING HAIR
- ✓ REVITALIZES DAMAGED HAIR
- ✓ EASY TO USE & NO KNOWN ADVERSE SIDE-EFFECTS

WWW.CAPILLUS.COM • TOLL-FREE USA: 1 (888) 272-9599



MOBILE LASER THERAPY CAP - BY PHYSICIANS ONLY

IS CAPILLUS LASER THERAPY RIGHT FOR YOU?

If your hair loss is due to heredity, Capillus laser therapy may be a great treatment option for you. Hair loss is a serious issue that affects millions of men and women around the world. The most common type of hair loss is a result of your genes: androgenetic alopecia (pattern baldness). It's a progressive condition (worsens over time) and constitutes 98% of hair loss. If you begin treatment during the earlier stages of progression, the CapillusRX™ may stop further progression and regrow your thinning hair.


HOW DOES LASER THERAPY REGROW HAIR?

Just as light makes plants grow, it helps hair regrow. Low-level laser therapy (LLLT) is based on the biological effects of photobiomodulation in living organisms and has been shown to be effective in stimulating and energizing the cells within the hair follicle. Hair Loss in androgenetic alopecia occurs through the process of miniaturization, where hair follicles that were formerly producing healthy hairs begin to produce thinner, shorter, more brittle hairs with weaker shafts until they eventually die and cease to produce hair permanently. Capillus laser therapy reverses the natural process of miniaturization on viable follicles, stimulating and energizing the cells within the hair follicle renewing the follicle for thicker, healthier hair.

CLINICALLY PROVEN TREATMENT

Capillus's independently-reviewed, double-blind clinical trial was registered on ClinicalTrials.gov and was managed by a neutral third party. The clinical trial data indicated that low-level laser therapy of the scalp every other day for 17 weeks with the CapillusPro™ device, significantly improved hair counts in those study participants who used the active (non-placebo) device. Subjects were able to use the device on a self-treatment home-use basis and no adverse events or side-effects were reported. The CapillusPro™ was proven to be a safe and effective treatment for androgenic alopecia. Results of the clinical study were recently published in a highly esteemed medical journal, *Dermatologic Surgery*. See details of clinical trials at: www.clinicaltrials.gov

PATIENT CASES



RESULTS MAY VARY



RESULTS MAY VARY

FEMALE BEFORE & AFTER

MALE BEFORE & AFTER

MOBILE LASER THERAPY - BY PHYSICIANS ONLY

Maximum Coverage | More Noticeable Results, Fastest

DISCREET THERAPY

You can wear your CapillusRX™ without embarrassment, any time, anywhere. Treat your hair loss discreetly under the included sports cap, or use under any hat of your choice.

EXCELLENT COVERAGE

Lasers work under the cap during your 6-minute treatment session, while you go about your everyday activities, providing light energy to the hair follicles at the cellular level. This cap features all laser technology (Class 3a).

NEW COMFORTABLE, FLEXIBLE FIT

A soft, flexible interior adapts to different head shapes and sizes, making Capillus laser therapy more comfortable than ever.



FAST TREATMENT, ANY TIME, ANYWHERE

Wear just 6 minutes every day. That's the time it takes to drink your morning coffee. You can use at home or anywhere else.



JUST 6 MINUTES A DAY



Laser Density

1560

mW Total Output

5

Year Warranty



WEAR IT ANYTIME, ANYWHERE.



Effective Treatment

“ It’s kind of like sending your hair to the gym. It’s a bio-stimulatory therapy, in that it stimulates the hair cells to both stay in the growth phase longer, reenter the growth phase sooner, and stay out of the non-growth phase for as long as possible. ”

Dr. Ken Anderson, MD, Atlanta, Hair Restoration Surgeon

Wear and Go

“ I like the type of laser device that you can actually put on your head and walk around because it has a battery pack. That seems to be the laser that is most compliant with patients and compliance means efficacy. ”

Dr. Shelly Friedman, DO, Scottsdale, Hair Restoration Surgeon

Scientifically Proven

“ What’s important to me about this treatment is that, 1) it is FDA cleared (that is of utmost importance), and 2) there are scientific studies that show its effectiveness treating a male and female pattern hair loss. ”

Dr. Bernard Nusbaum, MD, Miami, Hair Restoration Surgeon

¹ PHYSICIANS IN THE CAPILLUS NETWORK RECOMMEND CAPILLUS LASER THERAPY FOR THEIR PATIENTS.

² CAPILLUS DEVICES ARE MANUFACTURED IN THE USA BY HAND FROM FOREIGN SOURCED RAW MATERIALS.

³ LIMITED WARRANTY: WE PROVIDE A 5 YEAR LIMITED WARRANTY ON THE CAPILLUSRX™ LASER CAP DEVICE AND A 1 YEAR LIMITED WARRANTY ON THE BATTERY PACK & AC ADAPTER. CAPILLUS WARRANTS THE CAPILLUSRX™ AGAINST DEFECTS IN MATERIALS OR WORKMANSHIP FOR A PERIOD OF 5 YEAR FROM THE DATE OF PURCHASE FROM A CAPILLUS AUTHORIZED CENTER. TO ACTIVATE YOUR WARRANTY ON THE CAPILLUSRX™, REGISTER YOUR DEVICE WITHIN 30 DAYS OF PURCHASE AT WWW.CAPILLUS.COM/ACTIVATEWARRANTY. IF A DEFECT IN MATERIAL OR WORKMANSHIP EXISTS, CAPILLUS WILL REPAIR OR REPLACE THE NECESSARY PARTS WITH NEW OR REWORKED PARTS AT NO CHARGE. ANY REPLACEMENT PART WILL BE COVERED UNDER THE ORIGINAL 5 YEAR LIMITED WARRANTY. THIS WARRANTY IS NOT TRANSFERABLE & DOES NOT COVER DAMAGES CAUSED BY ACTS OF GOD, MISUSE, NEGLIGENCE, ACCIDENT, OR MODIFICATION OF ANY PART OF THE UNIT AS DETERMINED BY OUR RETURN DEPARTMENT. THIS WARRANTY DOES NOT COVER DAMAGE DUE TO IMPROPER OPERATION OR MAINTENANCE, CONNECTION TO IMPROPER VOLTAGE SUPPLY, OR ATTEMPTED REPAIR BY ANYONE OTHER THAN THE COMPANY. FAILURE TO USE AND MAINTAIN THE DEVICE IN ACCORDANCE WITH THE INSTRUCTIONS IN THE MANUAL WILL VOID THE PRODUCT WARRANTY.

⁴ CLINICAL TRIAL DATA AVAILABLE AT CLINICALTRIALS.GOV IDENTIFIER: NCT01967277 FOR RX: RESULTS FOR CAPILLUSPRO DEVICE.

↑ LASER APERTURE ↑ CAUTION	 LASER RADIATION - AVOID DIRECT EYE EXPOSURE 5mW Max CW Output @ 650nm Complies with FDA performance standards for laser products except for conformance with IEC 60825-1 Ed. 3 and IEC 60601-2-22 Ed. 3.1, as described in Laser Notice No. 56, dated May 8, 2019. CLASS 3R LASER PRODUCT
--------------------------------------	--